



The NMRPipe System for Multidimensional NMR Processing and Analysis



UNIVERSITY OF MARYLAND | NIST
INSTITUTE FOR BIOSCIENCE
& BIOTECHNOLOGY RESEARCH



frank.delaglio@nist.gov
Mar 28 2018

Key Steps of Spectral Processing

- Window Function and First Point Scaling
- Zero Fill
- Fourier Transform
- Phase Correction

Other Common Spectral Processing Tasks

- Solvent Subtraction
- Baseline Correction
- Linear Prediction
- Hilbert Transform and Inverse Processing

Special Cases





- Gradient-Enhanced Data
- Interleaved Experiments and Spectral Series
- Non-Uniform Sampling and Alternatives to Fourier Transform
- Correction of Bad Points

Post-Processing

- Generating Projections and Extracts from 3D or 4D Spectra
- Viewing 2D Strips from 3D Spectra
- Peak Detection and Fitting of Spectra and Spectral Series
- Extracting Chemical Shift Evolutions from Spectral Series
- Backbone Structure from Chemical Shifts and Dipolar Couplings
- Multivariate Analysis of Spectral Series
- Structure Manipulation and Analysis

With Thanks to Our NMRPipe Collaborators

frank.delaglio@nist.gov



National Institutes
of Health


Karen Allen • Luke Arbogast • Yves Aubin • Andrew Baldwin
Joseph Barchi • Ad Bax • Paul Bowyer • Robert Brinson
James Chou • Gabriel Cornilescu • Stefano Ciurli • Brian Coggins
Kathleen Farley • David Fushman • Michelle Gill • George Gray
Alex Grishaev • Stephan Grzesiek • Jeff Hoch • Sven Hyberts
Brad Jordan • Lewis Kay • Georg Kontaxis • Krish Krishnamurthy
John Kuszewski • Dong Long • Massimo Lucci • Mark Maciejewski
John Marino • Ryan McKay • Leszek Poppe • John Pfeiffer
Ben Ramirez • David Rovnyak • Dave Russell • Michael Shapiro
Adam Schuyler • Evgeny Tishchenko • Desi Tsao • Tobias Ulmer
Gerteen Vuister • Greg Walker • Gerard Weatherby • Justin Wu
Shen Yang • Jinfa Ying • Edward Zartler
Guang Zhu • Markus Zweckstetter



Frank Delaglio

NMRDraw Graphical Interface

Interactive Processing, Phasing and Peak Detection
1D-4D Data



New: NMRDraw Graphical Interface for Processing Script Generation

Created FT Script nmr_ft.com (on urz-testlab-02-1.public.unibas.ch)

File Draw 1D Mouse Peak Proc 2D File: ./test.fid (1603x2) (170x2)

X: 1 Y: 1 Z: 1 A: 1 Plus: 10 Minus: 10 First: 1156.31 Factor: 1.30

P0: P1: P0: -51.0 P1: 0.0 Pivot: 1 Phasing: On Off H-Axis: HN

Update Script Edit Script Execute Script

Show FID Show FT

X-Axis: HN
Solvent Filter: Convolution (S0L)

Sine Offset: 0.5 Endpoint: 0.95 Power: 2
Exp Hz: 0.0 Gauss Hz: 0.0 Gauss Offset: 0.0

Zero Fill: 1 Round to Power of 2:

P0: -51.0 P1: 0.0

Region Start: 3% End: 47%

Baseline Correction: Order 4

Y-Axis: 15N
Sine Offset: 0.5 Endpoint: 0.95 Power: 1
Exp Hz: 0.0 Gauss Hz: 0.0 Gauss Offset: 0.0

Zero Fill: 1 Round to Power of 2:


FT: Reverse (-neg) Rotate Halves (-alt)

P0: 0.0 P1: 0.0

Linear Prediction Extrapolation (LP):
Use Y-Axis Linear Prediction:
Y-Axis LP Order: 8 LP Mode: Forward-Backward

Special:
NUS-Style Extrapolation: None IST SMILE
Other Options:

FID Name: test.fid




Created FT Script nmr_ft.com

15N Pts Y


HN Pts X

Applications Places System [Terminal] Terminal delaglio_1@urz-testlab... Terminal Created FT Script nmr_... Fri Aug 4, 5:56 PM Frank Delaglio

Uniformly Sampled Data and 33% NUS Reconstructed by NMRPipe's Iterative Soft Thresholding (IST)



NUS Zero Fill as an Alternative to Linear Prediction for Conventional NMR Data (4D Methyl-Methyl NOE)



Key Steps of Spectral Processing


- Window function and first point scaling
- Zero Fill
- Fourier Transform
- Phase Correction

Spectral Processing Function as a UNIX Filter

```
nmrPipe -fn FT < test.fid > test.ft1
```

... or ...

```
nmrPipe -fn FT -in test.fid -out test.ft1
```




Spectral Processing Scheme as a UNIX Pipeline



```
nmrPipe -in test.fid -fn SP | nmrPipe -fn FT -out test.ft1
```

... Or ...




```
nmrPipe -in test.fid \  
| nmrPipe -fn SP \  
| nmrPipe -fn FT \  
-out test.ft1
```



Multidimensional Spectral Processing Schemes as a UNIX Pipelines



Related Software from the Bax Group at the NIH: TALOS-N Prediction of Protein Backbone and Sidechain Angles from Chemical Shifts








Yang Shen, and Ad Bax, *J. Biomol. NMR*, 56, 227-241(2013).

Related Software from the Bax Group at the NIH: SPARTA+ Prediction of Protein Chemical Shifts from Structure




SPARTA+ 2011205101229_001164_852904 <http://spin.niddk.nih.gov/bax/nmrserver/sparta>



Yang Shen and Ad Bax, *J. Biomol. NMR*, 48, 13-22 (2010)




2D HN/N NMR of NISTmAb as Simulated via *DYNAMO*, *SHIFTX2*, and *SPARTA+*




NMRPipe DC: Modules for Analysis of Dipolar Couplings






DC 2011205101202_001162_453935 <http://spin.niddk.nih.gov/bax/nmrserver/dc>



Tensor Parameters:
 Da: -6.8589e-04 Dr: -3.2939e-04 GMag: 2.1011e-03 Rh: 0.480
 Euler Rotations, Degrees: X: 18.53 Y: -27.60 Z: -55.64
 Order Matrix Values:
 -8.7887e-04 -1.0887e-04 -6.7871e-04 -7.4130e-04 -2.5099e-04
 sXX: 1.9180e-04 sYY: 1.1800e-03 sZZ: -1.3718e-03






DC 2011205101202_001162_453935 <http://spin.niddk.nih.gov/bax/nmrserver/dc>



Tensor Parameters:
 Da: -6.8589e-04 Dr: -3.2939e-04 GMag: 2.1011e-03 Rh: 0.480
 Euler Rotations, Degrees: X: 18.53 Y: -27.60 Z: -55.64
 Order Matrix Values:
 -8.7887e-04 -1.0887e-04 -6.7871e-04 -7.4130e-04 -2.5099e-04
 sXX: 1.9180e-04 sYY: 1.1800e-03 sZZ: -1.3718e-03

NMRPipe Protein Structure Analysis and Report Generation




PDB SS Utility 2012026140916_031907_502954 <http://spin.niddk.nih.gov/bax/nmrserver/pdbutil>





3M4Z.pdb Residues 1 to 309

Atoms: 2716 Mass: 36383 Da Area: 14921.1 A² Volume: 40847.2 A³ Radius of Gyration: 19.91 A



Distance Min: 1.20 A Max: 66.06 A Range X: -20.47 to 33.41 A Y: -35.39 to 22.92 A Z: -35.95 to 10.92 A Center XYZ: 7.16 -5.50 -13.60 A





NMRPipe Protein Structure Analysis and Report Generation



NMRPipe Spectrum versus Spectrum Correlation and Statistics





NMRPipe 1D Spectral Series Viewing with Principal Component Analysis (PCA)



Exploratory Analysis of 11 Varieties of Tea




NMRPipe 2D Spectral Series Viewing with Principal Component Analysis (PCA)



Monitoring Effects of Excipients, Formulation Parameters and Mutations on the High Order Structure of Filgrastim by NMR. Pharm Res, DOI 10.1007/s11095-015-1713-3 (2015)

Yves Aubin, Derek J. Hodgson, William B. Thach, Geneviève Gingras, and Simon Sauvé




NMRPipe Extraction of and Fitting of Chemical Shift Titration Curves






An NMRPipe Prototype Application: PCA Applied to Hyperspectral Optical Imaging

File Acquire Process


Help



0.0 Brightness Level Min
100.0 Brightness Level Max
125 Spectral Plane



Importing Spectral Data into a Virtual Reality Environment



About the Cover

November 4, 2010: Vol. 114, Iss. 43

Conventional (1D) electronic spectra of NO_2 are largely patternless, but coherent 2D spectra of NO_2 show numerous X-shaped clusters such as that shown here. This spectrum is displayed in 2D contour form (below, in purple and black) and as a corresponding 3D intensity surface (above, in purple and blue). Orthogonal 1D projections (orange) of the surface resemble conventional 1D spectra. The spectral data was preprocessed using a pipeline-based software system for multidimensional spectroscopy (Delaglio, F.; Grzesiek, S.; Vuister, G. W.; Zhu, G.; Pfeifer, J.; Bax, A. J. *Biomol. NMR* **1995**, *6*, 277–293), which was also used to generate images used to build the 3D scene. The figure was produced by Frank Delaglio (www.nmrscience.com), who used the interactive virtual reality environment Second Life to render this image (see page 11365).

Typical NMRPipe Workflow

The typical NMRPipe processing workflow involves frequent alternation between the UNIX command line and the `nmrDraw` graphical interface:

1. Prepare a script to convert spectrometer format data to NMRPipe format via an interactive interface. Commands: `bruker` / `varian` / `delta`
2. Execute the conversion script, and inspect the result via the `nmrDraw` graphical interface. Use the `Proc/Auto Process 1D` menu option to process the first 1D vector, and use the graphical interface to find a phase correction value.
3. At the UNIX command-line, use a text editor to create or adjust a processing script containing the needed phase corrections etc. Commonly, the processing script is copied from a previous case. Processing scripts can also be created via new general-purpose scripts: `basicFT2.com` / `basicFT3.com` / `basicFT4.com`.
4. Execute the processing script and read the processed result with `nmrDraw`. View vectors from any of the dimensions, and choose new phase correction values as needed.
4. Repeat steps 3 and 4 to adjust processing parameters as needed.
5. In the case of 3D data, create and inspect projections and view strips plots of the 3D data to confirm that the data is processed and phased properly. For example:

```
proj3D.tcl -in ft/test%03d.ft3 -abs  
peakHN.tcl -in ft/test%03d.ft3 -out hn.proj.tab -proj -hi Full  
scroll.tcl -in ft/test%03d.ft3 -tab hn.proj.tab -pair -hi 8%
```


Processing Functions for *nmrPipe* -fn ...

SOL	Time-domain Solvent Subtraction by Convolution and Subtraction
POLY -time	Time-domain Solvent Subtraction by Polynomial Fitting and Subtraction
SP	Window Function and First Point Scaling (-c 0.5 or -c 1.0)
ZF	Zero Fill (-auto rounds to power of two)
FT	Fourier Transform
PS	Phase Correction (-di deletes imaginaries after)
EXT	Extraction of a Chemical Shift Range
POLY -auto	Automated Baseline Correction by Polynomial Fitting and Subtraction
HT	Hilbert Transform (reconstruct previously deleted imaginary data)
TP	Exchange X-Axis and Y-Axis (also called YTP) for 2D, 3D, or 4D Data
ZTP	Exchange X-Axis and Z-Axis for 3D and 4D Data
ATP	Exchange X-Axis and A-Axis for 4D Data
LP	Linear Prediction (extend or replace time-domain data)
MEM	Maximum Entropy Reconstruction
ML	Maximum Likelihood Frequency Map
SMILE	SMILE Reconstruction
COADD	Weighted Co-Addition of Adjacent Vectors for Interleaved Data
MAC	NMRPipe's Macro Language for Simple Custom Processing Functions

Listing Command Line Options: nmrPipe -help -fn ...

```
% nmrPipe -help -fn ZF
```

ZF: Extend By Zero Filling;

Use only one of the following:

-zf **zfCnt** **[1]** Number of Times to Double the Size.
-pad **padCnt** Zeros to Add.
-size **finSize** Final Size.

Other Flags:

-mid Zero Fill in Middle.
-inter Zero Fill by zfCount Between Points.
-auto Round Final Size to Power of 2.

Removing Previous Zero Filling:

-inv Extract Original Time Domain.

Set Time-Domain Size for Zero Fill Extrapolation:

-td Set Time-Domain Size to Zero Fill Size.
-notd Do Not Change Time-Domain Size (Default).

Typical NMRPipe 3D Spectral Processing Scheme

Commonly Adjusted Parameters / X-Axis Processing / Y-Axis Processing / Z-Axis Processing

The Script Provides Reproducible Processing

```
#!/bin/csh
```

```
xyz2pipe -in fid/test%03d.fid -x -verb \ Read vectors from X-Axis
| nmrPipe -fn SOL \ Solvent Subtraction
| nmrPipe -fn SP -off 0.5 -end 0.98 -pow 2 -c 0.5 \ Window and First Point Scale
| nmrPipe -fn ZF \ Zero Fill
| nmrPipe -fn FT \ Fourier Transform
| nmrPipe -fn PS -p0 43.0 -p1 0.0 -di \ Phase Correction
| nmrPipe -fn EXT -x1 10.5ppm -xn 5.7ppm -sw \ Extract PPM Range
| nmrPipe -fn TP \ X/Y Transpose
| nmrPipe -fn SP -off 0.5 -end 0.98 -pow 1 -c 1.0 \ Window and First Point Scale
| nmrPipe -fn ZF \ Zero Fill
| nmrPipe -fn FT \ Fourier Transform
| nmrPipe -fn PS -p0 -90.0 -p1 180.0 -di \ Phase Correction
| nmrPipe -fn TP \ X/Y Transpose
| nmrPipe -fn POLY -auto \ Baseline Correction
| pipe2xyz -out ft/test%03d.ft2 -x Write vectors to X-Axis

xyz2pipe -in ft/test%03d.ft2 -z -verb \ Read vectors from Z-Axis
| nmrPipe -fn SP -off 0.5 -end 0.98 -pow 1 -c 0.5 \ Window and First Point Scale
| nmrPipe -fn ZF \ Zero Fill
| nmrPipe -fn FT \ Fourier Transform
| nmrPipe -fn PS -p0 0.0 -p1 0.0 -di \ Phase Correction
| pipe2xyz -out ft/test%03d.ft3 -z Write vectors to Z-Axis
```

New General-Purpose Processing and Script Generation:

`basicFT2.com basicFT3.com basicFT4.com`

The latest version of NMRPipe includes several new general-purpose commands for simplified spectral processing:

- Default input file is based on the contents of the current directory, so you can get a quick and easy trial result just with a command like:

```
basicFT2.com -xP0 43
```

- Default parameters change automatically for amide-detected data.
- There are options for **Linear Prediction** and **NUS Zero Fill** (IST reconstruction applied to conventional data, described later).
- First-point time-domain scaling is adjusted automatically according to first-order phase correction.
- The commands can perform processing directly, or **generate the NMRPipe script** so that you can edit it yourself.
- Command `basicFT2.com` can process a **2D spectrum**, extract and process a **2D XY or XZ plane** from 3D or 4D data, a **2D XA plane** from 4D data, or a process an entire pseudo-3D series of 2D spectra.
- `basicFT3.com` and `basicFT4.com` process **3D and 4D** spectra.

New General-Purpose Processing and Script Generation: *basicFT2.com basicFT3.com basicFT4.com*

This:

```
basicFT2.com -xP0 13 -xP1 0 -xEXTX1 10.4ppm -xEXTXN 5.4ppm -yP0 -90 -yP1 180
```

*Input (-in) determined
automatically*

takes the place of this:


```
nmrPipe -in fid/test001.fid \  
| nmrPipe -fn SOL \  
| nmrPipe -fn SP -off 0.5 -end 0.98 -pow 2 -c 0.5 \  
| nmrPipe -fn ZF -auto \  
| nmrPipe -fn FT \  
| nmrPipe -fn PS -p0 13 -p1 0.0 -di \  
| nmrPipe -fn EXT -x1 10.4ppm -xn 5.4ppm -sw \  
| nmrPipe -fn TP \  
| nmrPipe -fn SP -off 0.5 -end 0.98 -pow 1 -c 1.0 \  
| nmrPipe -fn ZF -auto \  
| nmrPipe -fn FT \  
| nmrPipe -fn PS -p0 -90 -p1 180 -di \  
| nmrPipe -fn TP \  
| nmrPipe -fn POLY -auto \  
-out test.ft2 -verb -ov
```

SOL default for HN data


*First-Point Scale Factors
(-c values) automatically
calculated from P1*

*POLY -auto default for
HN data*

Output (-out) set by default



Non-Uniform Sampling in NMRPipe



About Non-Uniform Sampling

Non-Uniform Sampling (NUS) is an acquisition method for multidimensional NMR spectra that works by skipping some fraction of the data that would be acquired in a conventional measurement, which is uniformly sampled (US).

The goal of NUS is to **improve the spectral resolution** (or other metric of spectral quality) obtained with a given amount of measurement time.

Because some of the data is skipped in a NUS acquisition, the usual Fourier transform processing used for conventional data is not ideal, and so other reconstruction methods are required to take best advantage of NUS.

There are many tools for reconstructing NUS data, including Matrix Decomposition (MDD software, Orekhov), Maximum Entropy (RNMRTK, Hoch et al.), Iterative Soft Thresholding (IST) (hmsIST, Hyberts and Wagner, and applications from Stern, Donoho, Hoch et al.), SCRUB (Coggins and Zhou), and NESTA (Sun, Gill, Byrd et al.). Many of these tools use NMRPipe as part of the reconstruction workflow.

NMRPipe includes its own implementation of IST, so that it is easy to perform all the steps needed for NUS reconstruction. This presentation explains NUS and IST, and shows how to **use NMRPipe to convert and process 2D and 3D NUS** NMR data with IST. It also shows how to apply IST to conventional data as an **alternative to Linear Prediction**.

NUS Software from Our Friends

Many NUS Reconstruction Systems Use NMRPipe Workflows



SMILE spin.niddk.nih.gov/bax/software/smile



NESTA nestanmr.com




SCRUB coggins.biochem.duke.edu/scrub



hmsIST gwagner.med.harvard.edu/intranet/istHMS

*Find these methods and
more on NMRbox*



NMRbox.org

Non-Uniform Sampling in NMRPipe

The NMRPipe NUS tools expect the following:

- NUS schedules are plain text files with the same number of space-separated values on each line, and with no blank lines or other content. Increment numbers are always given as integers.
- 2D NUS schedules have one increment number per line.
- 3D NUS schedules have two increment numbers per line, one for the Y-Axis, and one for the Z-Axis.
- 4D schedules have three increment numbers per line, for the Y-Axis, Z-Axis, and A-Axis.

There is more than one convention for recording NUS sampling schedules. NMRPipe can directly accommodate the following:

- Some schedules will use increment numbers which start at one for the first increment rather than zero.
- Some schedules will include floating-point weighting factors along with the integer increment numbers.
- By default, NMRPipe assumes that 3D and 4D schedules have the Y-Axis increment value listed first on each line. In some conventions, this order is reversed, so that the Y-Axis increment is last on each line.

NMRPipe Workflow for NUS Data

1. Prepare a conversion script with the usual tools, including options which sort and expand the NUS time-domain data with zeros. Commands:

```
bruker -nus / varian -nus / delta -nus
```

2. Perform **trial processing** via ordinary Fourier processing on the NUS fid, to confirm phase correction values, etc. This can be done via the usual NMRPipe scripts, or via new general-purpose processing scripts:

```
basicFT2.com / basicFT3.com / basicFT4.com
```

3. In the case of 3D or 4D data, inspect the ordinary Fourier transform result by preparing projections. Commands:

```
proj3D.tcl / proj4D.tcl
```

4. Use the ordinary Fourier transform result to estimate the convergence parameter used by NMRPipe's NUS reconstruction tools:

```
specStat.com -stat istMaxRes -in ...
```

5. Use NMRPipe's NUS reconstruction scripts to generate a final spectrum:

```
ist2D.com / ist3D.com / ist4D.com
```

How to Convert NUS Data in NMRPipe

NMRPipe Conversion Utility Version 2015.244.10.57

Spectrometer Input: /ser
NUS Schedule: /nuslist
Output Template: /fid/test%03d.fid
NUS Mask Output: /mask/test%03d.fid
Output Script: fid.com
Other Options:

Input Protocol: Bruker NUS
Output Protocol: NMRPipe
Dimension Count: N
2D Mode: States
Temperature (K): From File
NUS Samples: Auto
NUS Index Offsets: Auto
 Reverse NUS Column Order

Digital Oversampling Correction: During Conversion (Normal FID) During Processing (Better Baseline)

	x-axis	y-axis	z-axis
Total Points R+I:	1024	128	64
Valid Points:	512	64	32
Acquisition Mode:	Complex	Complex	Complex
Spectral Width Hz:	10000.00	10000.00	2000.000
Observe Freq MHz:	500.000	500.000	100.00
Center Position PPM:	H2O	4.700	50.00
Axis Label:	X	Y	Z

Read Parameters Save Script Execute Script Hide Script Clear Script Update Script Quit Help

1. Select the binary spectrometer-format input data as usual.

2. Select the NUS Sampling Schedule Input File.

3. Choose options for best baseline.

4. Read Parameters
This will also read the NUS schedule, and extract information about sample count, dimension sizes, etc.

Bruker 3D NUS Conversion

```
#!/bin/csh
```

```
nusExpand.tcl -mode bruker -sampleCount 2048 -off 0 \  
-in ./ser -out ./ser_full -sample ./nuslist
```

```
bruk2pipe -in ./ser_full \  
-bad 0.0 -aswap -AMX -decim 1680 -dspfvs 20 -grpdlly 67.9866027832031 \  
-xN 2048 -yN 64 -zN 32 \  
-xT 1024 -yT 32 -zT 160 \  
-xMODE DQD -yMODE Echo-AntiEcho -zMODE Complex \  
-xSW 11904.762 -ySW 3846.154 -zSW 11904.762 \  
-xOBS 950.204 -yOBS 96.294 -zOBS 950.204 \  
-xCAR 4.773 -yCAR 118.579 -zCAR 4.773 \  
-xLAB HN -yLAB 15N -zLAB 1H \  
-ndim 3 -aq2D States \  
-out ./fid/test%03d.fid -verb -ov
```

```
xyz2pipe -in ./fid/test%03d.fid \  
| nusExpand.tcl -mask -noexpand -mode pipe -sampleCount 2048 -off 0 \  
-in stdin -out ./mask/test%03d.fid -sample ./nuslist
```

Expand the Spectrometer-format Data / Convert Expanded Data / Create Mask

Varian 3D NUS Conversion

```
#!/bin/csh
```

```
nusExpand.tcl -mode varian -sampleCount 455 -off 0 \  
-in ./fid -out ./fid_full -sample ./sampling.sch
```

```
var2pipe -in ./fid_full \  
-noaswap -aqORD 1 \  
-xN 2404 -yN 140 -zN 52 \  
-xT 1202 -yT 70 -zT 26 \  
-xMODE Complex -yMODE Complex -zMODE Rance-Kay \  
-xSW 12019.231 -ySW 7000.000 -zSW 1320.000 \  
-xOBS 599.422 -yOBS 599.422 -zOBS 60.746 \  
-xCAR 4.770 -yCAR 4.770 -zCAR 118.178 \  
-xLAB HN -yLAB 1H -zLAB N15 \  
-ndim 3 -aq2D States \  
-out ./data/test%03d.fid -verb -ov
```

```
xyz2pipe -in ./data/test%03d.fid \  
| nusExpand.tcl -mask -noexpand -mode pipe -sampleCount 455 -off 0 \  
-in stdin -out ./mask/test%03d.fid -sample ./sampling.sch
```

Expand the Spectrometer-format Data / Convert Expanded Data / Create Mask

JEOL 3D NUS Conversion

```
#!/bin/csh
```


```
delta2pipe -in ./3C15N_Ubq_hnco_20NUS-1-1.jdf -nusDim 3 \  
-nodf -dfVal 19.687500 -trVal 8.000000e-01 \  
-xN 2048 -yN 1004 -zN 1 \  
-xT 1024 -yT 50 -zT 25 \  
-xMODE Complex -yMODE Complex -zMODE Real \  
-xSW 11261.261 -ySW 3018.594 -zSW 3041.363 \  
-xORIG -2815.707 -yORIG 25428.789 -zORIG 5783.051 \  
-xOBS 600.172 -yOBS 150.913 -zOBS 60.815 \  
-xCAR 4.681 -yCAR 176.000 -zCAR 120.000 \  
-xFT Time -yFT Time -zFT Time \  
-xLAB HN -yLAB CO -zLAB 15N \  
-ndim 2 -aq2D Complex \  
-out nus.fid -verb -ov
```

```
nusExpand.tcl -in nus.fid -out fid/test%03d.fid \  
-sample hnco_NUS.txt -sampleCount Auto
```


```
nusExpand.tcl -in nus.fid -out mask/test%03d.fid \  
-sample hnco_NUS.txt -sampleCount Auto -mask
```

Convert 3D Data as NUS 2D / Expand the Converted Data to 3D / Create Mask

Output of NMRPipe NUS Conversion



*Sorted, expanded time-domain data,
arranged like conventional data but
with the increments that were skipped
filled in with zeros*



*Corresponding mask of ones and
zeros, used as additional input for
NMRPipe reconstruction tools*

How to Process NUS Data in NMRPipe

It is recommended to **apply trial processing** of the expanded time-domain NUS data using ordinary Fourier transform schemes. This often allows confirmation of details such as phase correction before more time-consuming reconstruction methods are used. General-purpose scripts basicFT2.com basicFT3.com basicFT4.com are provided to help with this step ...

```
#!/bin/csh
```

```
basicFT3.com \  
-in fid/test%03d.fid -out ft/test%03d.ft3 \  
-xP0 -75 -xP1 0 -xEXTX1 10.4ppm -xEXTXN 5.4ppm \  
-zFTARG alt
```

Process complete 3D

```
proj3D.tcl -in ft/test%03d.ft3 -abs
```

*Generate and Inspect 2D
projections from 3D*

How to Process NUS Data with NMRPipe's Iterative Soft Thresholding (IST)

After processing parameters have been confirmed using ordinary Fourier transform schemes, use the NMRPipe IST script ([ist2D.com](#) [ist3D.com](#) [ist4D.com](#)) to reconstruct a spectrum.



For convenience, the NMRPipe IST scripts have similar arguments to scripts [basicFT2.com](#) etc. NMRPipe's IST takes the expanded **time-domain data** and corresponding **mask** as input, and produces a **spectrum** as output.

The parameter [istMaxRes](#) defines the convergence condition for IST, specified as a percentage of the largest value in the starting spectrum. Convergence is achieved when the largest value in the residual is below this percentage. By default, this parameter is determined automatically, and all the examples in the NMRPipe demo data archive use the automated setting.



```
ist3D.com -istMaxRes Auto \  
-in fid/test%03d.fid -mask mask/test%03d.fid -out ist/test%03d.ft3 \  
-xP0 -75 -xP1 0 -xEXTX1 10.4ppm -xEXTXN 5.4ppm \  
-zFTARG alt
```


NUS Zero Fill as an Alternative to Linear Prediction

Uniformly Sampled Data



Non-Uniformly Sampled Data



- Conventional Fourier transform, Conventional Uniform Sampling Schedule: Broad lines and periodic truncation wiggle artifacts
- Conventional Fourier transform, Non Uniform Sampling Schedule: Narrow lines and random-noise-like artifacts

All peaks for a given schedule have the same artifacts, and the size of the artifacts is proportional to the size of the peaks.


As noted by Hoch and coworkers, IST can be used with Non-Uniform schedules, and also with conventional Uniform Schedules. This means it is possible to apply IST to conventionally-sampled data as an alternative to extrapolating data via Linear Prediction. We call this approach **NUS Zero Fill**.

NUS Zero Fill using IST as an Alternative to Linear Prediction for Conventional NMR Data (^{15}N -NOE)


Conventional FT
4 sec

Linear Prediction
2 min 30 sec

NUS Zero Fill via IST
20 min 21 sec



NMRPipe is Provided by the National Institute of Standards and Technology in Partnership with the National Institutes of Health



NIST
**National Institute of
Standards and Technology**
U.S. Department of Commerce




UNIVERSITY OF MARYLAND | NIST
**INSTITUTE FOR BIOSCIENCE
& BIOTECHNOLOGY RESEARCH**




**National Institutes
of Health**

NMRPipe: a multidimensional spectral processing system based on UNIX pipes. Frank Delaglio, Stephan Grzesiek, Geerten. W. Vuister, Guang Zhu, John Pfeifer, and Ad Bax, *J. Biomol. NMR.* **6**, 277-293 (1995).




UMD | NIST
INSTITUTE FOR
BIOSCIENCE &
BIOTECHNOLOGY
RESEARCH

www.ibbr.umd.edu



www.nmrpipe.com



NIST Disclaimer: Certain commercial equipment, instruments, and materials are identified in this presentation in order to specify the experimental procedure. Such identification does not imply recommendation or endorsement by the National Institute of Standards and Technology, nor does it imply that the material or equipment identified is necessarily the best available for the purpose.